

MEIKO
Maschinenbau GmbH & Co. KG
Englerstraße 3

D- 77652 Offenburg

Phone: +49 (0) 385 5682 65
Fax: +49 (0) 385 5682 67
Email: info@hygcn.de
Web: www.hygcn.de

29.03.2018

Expertise

As the reports SN 8009 of 2008-07-01 and SN 8009a of 2008-07-21 show the bedpan washer TopLine 20 has been tested regarding the efficacy of disinfection. The test soiling RAMS (bovine albumin- mucin- maize starch) or dirty conditions (3g BSA/l, 3ml erythrocyte/liter) were mixed with spores of *Clostridium difficile* ATCC 6989. After exposition in the process in TopLine 20 no growth was detectable, and therefor a sufficient sporicidal disinfection in correlation to the A_0 -value was proved.

The further technical development of the badpan washer to TOPIC shows no changes affecting the efficacy of disinfection, why the results of the tests (SN 8009 and SN 8009a) apply also for the machine TOPIC.



Prof. Dr. med. H.-P. Werner